

TSI Model 3332 Dilution System

Rental and Application Notes

- Dilution ratios of 100 to 1 or 10 to 1
- Higher ratios achieved by setting up two diluters in tandem
- Very low particle loss from 0.5 to 10 μm
- 1 L flow rate
- The main application of the Model 3332 Dilution System is in the dilution of test aerosols for the Model 3330 Optical Particle Sizer.
- Other applications having a 1 liter per minute flow rate can also use the Model 3332 Particle Diluter



	Specifications
Weight	<5lbs
Dimensions	160 mm x 300 mm x 110 mm (6.3 in. x 11.8 in. x 4.3 in)
Power	9-15 VDC (AC adapter, battery operation)
Flow	1Lpm